

Product Datasheet

Mouse MBP ELISA Kit (Ready to Use) (orb553713)

Description	Mouse Major Basic Protein (MBP) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for c...
Reactivity	Mouse
Range	0.625-40ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q61878
Sensitivity	0.235ng/mL
Entrez	19074
Expiration Date	Please enquire.