

Product Datasheet

Mouse MASP2 ELISA Kit (Ready to Use) (orb553715)

Description	Mouse Mannose Associated Serine Protease 2 (MASP2) Sandwich ELISA KitReady to Use: Pre-diluted detec...
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Sensitivity	0.067ng/mL
Expiration Date	Please enquire.