

## Product Datasheet

### Mouse MCL1 ELISA Kit (Ready to Use) (orb553775)

|                          |                                                                                                                    |
|--------------------------|--------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>       | Mouse Myeloid Cell Leukemia Sequence 1, Bcl2 Related (MCL1) Sandwich ELISA KitReady to Use: Pre-dilu...            |
| <b>Reactivity</b>        | Mouse                                                                                                              |
| <b>Range</b>             | 0.156-10ng/mL                                                                                                      |
| <b>Form/Appearance</b>   | Ready to Use                                                                                                       |
| <b>Storage</b>           | All reagents should be stored at -20°C upon receipt                                                                |
| <b>Note</b>              | For research use only                                                                                              |
| <b>Application notes</b> | Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids. |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates and other biological fluids.                                                     |
| <b>Uniprot ID</b>        | <b>P97287</b>                                                                                                      |
| <b>Sensitivity</b>       | 0.056ng/mL                                                                                                         |
| <b>Entrez</b>            | <b>17210</b>                                                                                                       |
| <b>Expiration Date</b>   | Please enquire.                                                                                                    |