

Product Datasheet

Mouse NPNT ELISA Kit (Ready to Use) (orb553805)

Description	Mouse Nephronectin (NPNT) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for conveni...
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q91V88
Sensitivity	0.051ng/mL
Entrez	114249
Expiration Date	Please enquire.