

## Product Datasheet

# Mouse Glucocorticoid Receptor ELISA Kit (Ready to Use) (orb553850)

|                          |                                                                                                                   |
|--------------------------|-------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>       | Mouse Nuclear Receptor Subfamily 3, Group C, Member 1 (Glucocorticoid Receptor) Sandwich ELISA KitRe...           |
| <b>Reactivity</b>        | Mouse                                                                                                             |
| <b>Range</b>             | 78.1-5000pg/mL                                                                                                    |
| <b>Form/Appearance</b>   | Ready to Use                                                                                                      |
| <b>Storage</b>           | All reagents should be stored at -20°C upon receipt                                                               |
| <b>Note</b>              | For research use only                                                                                             |
| <b>Application notes</b> | Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids. |
| <b>Sample Types</b>      | tissue homogenates, cell lysates and other biological fluids.                                                     |
| <b>Sensitivity</b>       | 33.7pg/mL                                                                                                         |
| <b>Expiration Date</b>   | Please enquire.                                                                                                   |