

Product Datasheet

Mouse PSGL1 ELISA Kit (Ready to Use) (orb553873)

Description	Mouse P-Selectin Glycoprotein Ligand 1 (PSGL1) Sandwich ELISA KitReady to Use: Pre-diluted detection...
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q62170
Sensitivity	0.065ng/mL
Entrez	20345
Expiration Date	Please enquire.