

Product Datasheet

Mouse PRDX6 ELISA Kit (Ready to Use) (orb553896)

Description	Mouse Peroxiredoxin 6 (PRDX6) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for con...
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	008709
Sensitivity	0.063ng/mL
Entrez	11758
Expiration Date	Please enquire.