

Product Datasheet

Mouse PCK1 ELISA Kit (Ready to Use) (orb553911)

Description	Mouse Phosphoenolpyruvate Carboxykinase 1, Soluble (PCK1) Sandwich ELISA KitReady to Use: Pre-dilute...
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q9Z2V4
Sensitivity	0.058ng/mL
Entrez	18534
Expiration Date	Please enquire.