

Product Datasheet

Mouse uPAR ELISA Kit (Ready to Use) (orb553925)

Description	Mouse Plasminogen Activator, Urokinase Receptor (uPAR) Sandwich ELISA KitReady to Use: Pre-diluted d...
Reactivity	Mouse
Range	125-8000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P35456
Sensitivity	47pg/mL
Entrez	18793
Expiration Date	Please enquire.