

Product Datasheet

Mouse TPSb2 ELISA Kit (orb553938)

Description	Mouse Polyclonal Antibody to Tryptase Beta 2 (TPSb2) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Sensitivity	0.054ng/mL
Expiration Date	Please enquire.