

Product Datasheet

Mouse P4Ha1 ELISA Kit (orb553961)

Description	Mouse Prolyl-4-Hydroxylase Alpha Polypeptide I (P4Ha1) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Sensitivity	0.108ng/mL
Expiration Date	Please enquire.