

Product Datasheet

Mouse S100A4 ELISA Kit (Ready to Use) (orb554034)

Description	Mouse S100 Calcium Binding Protein A4 (S100A4) Sandwich ELISA KitReady to Use: Pre-diluted detection...
Reactivity	Mouse
Range	31.2-2000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P07091
Sensitivity	14.2pg/mL
Entrez	20198
Expiration Date	Please enquire.