

Product Datasheet

Mouse SCG2 ELISA Kit (Ready to Use) (orb554049)

Description	Mouse Secretogranin II (SCG2) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for con...
Reactivity	Mouse
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q03517
Sensitivity	0.113ng/mL
Entrez	20254
Expiration Date	Please enquire.