

Product Datasheet

Mouse S1PR1 ELISA Kit (Ready to Use) (orb554104)

Description	Mouse Sphingosine 1 Phosphate Receptor 1 (S1PR1) Sandwich ELISA KitReady to Use: Pre-diluted detecti...
Reactivity	Mouse
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	008530
Sensitivity	0.115ng/mL
Entrez	13609
Expiration Date	Please enquire.