

Product Datasheet

Mouse SGPP2 ELISA Kit (orb554105)

Description	Mouse Sphingosine-1-Phosphate Phosphatase 2 (SGPP2) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Sensitivity	0.073ng/mL
Expiration Date	Please enquire.