

Product Datasheet

Mouse SF3B3 ELISA Kit (Ready to Use) (orb554109)

Description	Mouse Splicing Factor 3B Subunit 3 (SF3B3) Sandwich ELISA KitReady to Use: Pre-diluted detection rea...
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q921M3
Sensitivity	0.056ng/mL
Entrez	101943
Expiration Date	Please enquire.