

Product Datasheet

Mouse STC1 ELISA Kit (Ready to Use) (orb554112)

Description	Mouse Stanniocalcin 1 (STC1) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for conv...
Reactivity	Mouse
Range	78.1-5000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	O55183
Sensitivity	29.3pg/mL
Entrez	20855
Expiration Date	Please enquire.