

Product Datasheet

Mouse SYNPO2 ELISA Kit (Ready to Use) (orb554137)

Description	Mouse Synaptopodin 2 (SYNPO2) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for con...
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	E9Q1U2
Sensitivity	0.051ng/mL
Entrez	118449
Expiration Date	Please enquire.