

Product Datasheet

Mouse THBS4 ELISA Kit (Ready to Use) (orb554160)

Description	Mouse Thrombospondin 4 (THBS4) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for co...
Reactivity	Mouse
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityplasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q9Z1T2
Sensitivity	0.125ng/mL
Entrez	21828
Expiration Date	Please enquire.