

Product Datasheet

Mouse VCC1 ELISA Kit (Ready to Use) (orb554265)

Description	Mouse VEGF Co Regulated Chemokine 1 (VCC1) Sandwich ELISA KitReady to Use: Pre-diluted detection rea...
Reactivity	Mouse
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q5UW37
Sensitivity	5.8pg/mL
Entrez	232983
Expiration Date	Please enquire.