

Product Datasheet

Mouse VGF ELISA Kit (Ready to Use) (orb554270)

Description	Mouse VGF Nerve Growth Factor Inducible (VGF) Sandwich ELISA KitReady to Use: Pre-diluted detection ...
Reactivity	Mouse
Range	31.2-2000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q0VGU4
Sensitivity	16.3pg/mL
Entrez	381677
Expiration Date	Please enquire.