

Product Datasheet

Mouse WASF2 ELISA Kit (Ready to Use) (orb554286)

Description	Mouse Wiskott Aldrich Syndrome Protein Family, Member 2 (WASF2) Sandwich ELISA KitReady to Use: Pre-...
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q8BH43
Sensitivity	0.053ng/mL
Entrez	242687
Expiration Date	Please enquire.