

Product Datasheet

Porcine ACE ELISA Kit (Ready to Use) (orb554306)

Description	Porcine Angiotensin I Converting Enzyme (ACE) Sandwich ELISA KitReady to Use: Pre-diluted detection ...
Reactivity	Porcine
Range	31.2-2000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	A1YQ85
Sensitivity	10.7pg/mL
Entrez	613133
Expiration Date	Please enquire.