

Product Datasheet

Porcine Apolipoprotein A1 ELISA Kit (Ready to Use) (orb554308)

Description	Porcine Apolipoprotein A1 (Apolipoprotein A1) Sandwich ELISA KitReady to Use: Pre-diluted detection ...
Reactivity	Porcine
Range	28.1-1800ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at - 20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P18648
Sensitivity	12.7ng/mL
Entrez	397691
Expiration Date	Please enquire.