

Product Datasheet

Porcine Apolipoprotein E ELISA Kit (Ready to Use) (orb554311)

Description	Porcine Apolipoprotein E (Apolipoprotein E) Sandwich ELISA KitReady to Use: Pre-diluted detection re...
Reactivity	Porcine
Range	1.25-80ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at - 20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	F1RM45
Sensitivity	0.56ng/mL
Entrez	397576
Expiration Date	Please enquire.