

Product Datasheet

Porcine BMP15 ELISA Kit (orb554314)

Description	Porcine Bone Morphogenetic Protein 15 (BMP15) Sandwich ELISA Kit
Reactivity	Porcine
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Sensitivity	10.4pg/mL
Expiration Date	Please enquire.