

Product Datasheet

Porcine FGF21 ELISA Kit (Ready to Use) (orb554349)

Description	Porcine Fibroblast Growth Factor 21 (FGF21) Sandwich ELISA KitReady to Use: Pre-diluted detection re...
Reactivity	Porcine
Range	31.2-2000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	C6JWG4
Sensitivity	13.3pg/mL
Entrez	100302504
Expiration Date	Please enquire.