

Product Datasheet

Porcine FGF7 ELISA Kit (Ready to Use) (orb554350)

Description	Porcine Fibroblast Growth Factor 7 (FGF7) Sandwich ELISA KitReady to Use: Pre-diluted detection reagent...
Reactivity	Porcine
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q9N198
Sensitivity	5.9pg/mL
Entrez	397281
Expiration Date	Please enquire.