

Product Datasheet

Porcine GSTp ELISA Kit (Ready to Use) (orb554358)

Description	Porcine Glutathione S Transferase Pi (GSTp) Sandwich ELISA KitReady to Use: Pre-diluted detection re...
Reactivity	Porcine
Range	62.5-4000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P80031
Sensitivity	22.8pg/mL
Entrez	100739163
Expiration Date	Please enquire.