

Product Datasheet

Porcine IGFBP1 ELISA Kit (Ready to Use) (orb554378)

Description	Porcine Insulin Like Growth Factor Binding Protein 1 (IGFBP1) Sandwich ELISA KitReady to Use: Pre-di...
Reactivity	Porcine
Range	3.12-200ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	E0A8P1
Sensitivity	1.32ng/mL
Entrez	397270
Expiration Date	Please enquire.