

Product Datasheet

Porcine IGFBP4 ELISA Kit (Ready to Use) (orb554381)

Description	Porcine Insulin Like Growth Factor Binding Protein 4 (IGFBP4) Sandwich ELISA KitReady to Use: Pre-di...
Reactivity	Porcine
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P24854
Sensitivity	0.098ng/mL
Entrez	100144490
Expiration Date	Please enquire.