

Product Datasheet

Porcine IL5 ELISA Kit (Ready to Use) (orb554395)

Description	Porcine Interleukin 5 (IL5) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for conve...
Reactivity	Porcine
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q9MYM5
Sensitivity	6.5pg/mL
Entrez	397409
Expiration Date	Please enquire.