

Product Datasheet

Porcine Leptin Receptor ELISA Kit (Ready to Use) (orb554400)

Description	Porcine Leptin Receptor (Leptin Receptor) Sandwich ELISA KitReady to Use: Pre-diluted detection reagent...
Reactivity	Porcine
Range	1.56-100ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	002671
Sensitivity	0.59ng/mL
Entrez	396836
Expiration Date	Please enquire.