

Product Datasheet

Porcine S100A12 ELISA Kit (Ready to Use) (orb554435)

Description	Porcine S100 Calcium Binding Protein A12 (S100A12) Sandwich ELISA KitReady to Use: Pre-diluted detec...
Reactivity	Porcine
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P80310
Sensitivity	0.043ng/mL
Entrez	100301483
Expiration Date	Please enquire.