

Product Datasheet

Porcine VEGFR1 ELISA Kit (Ready to Use) (orb554462)

Description	Porcine Vascular Endothelial Growth Factor Receptor 1 (VEGFR1) Sandwich ELISA KitReady to Use: Pre-d...
Reactivity	Porcine
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	F1RSU5
Sensitivity	0.115ng/mL
Entrez	396763
Expiration Date	Please enquire.