

Product Datasheet

Rabbit AGC ELISA Kit (Ready to Use) (orb554466)

Description	Rabbit Aggrecan (AGC) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for convenient ...
Reactivity	Rabbit
Range	3.12-200ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	G1U677
Sensitivity	1.12ng/mL
Entrez	100009079
Expiration Date	Please enquire.