

Product Datasheet

Rabbit Apolipoprotein AI ELISA Kit (Ready to Use) (orb554470)

Description	Rabbit Apolipoprotein A1 (Apolipoprotein AI) Sandwich ELISA KitReady to Use: Pre-diluted detection r...
Reactivity	Rabbit
Range	31.2-2000ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at - 20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P09809
Sensitivity	9.35ng/mL
Entrez	100009253
Expiration Date	Please enquire.