

Product Datasheet

Rabbit CNTF ELISA Kit (Ready to Use) (orb554480)

Description	Rabbit Ciliary Neurotrophic Factor (CNTF) Sandwich ELISA KitReady to Use: Pre-diluted detection reagent...
Reactivity	Rabbit
Range	31.2-2000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P14188
Sensitivity	11.8pg/mL
Entrez	100009117
Expiration Date	Please enquire.