

Product Datasheet

Rabbit IL1RA ELISA Kit (Ready to Use) (orb554498)

Description	Rabbit Interleukin 1 Receptor Antagonist (IL1RA) Sandwich ELISA KitReady to Use: Pre-diluted detecti...
Reactivity	Rabbit
Range	31.2-2000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P26890
Sensitivity	12.8pg/mL
Entrez	100009558
Expiration Date	Please enquire.