

Product Datasheet

Rabbit MIP1b ELISA Kit (Ready to Use) (orb554506)

Description	Rabbit Macrophage Inflammatory Protein 1 Beta (MIP1b) Sandwich ELISA KitReady to Use: Pre-diluted de...
Reactivity	Rabbit
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P46632
Sensitivity	6.2pg/mL
Entrez	100008984
Expiration Date	Please enquire.