

Product Datasheet

Rabbit MMP3 ELISA Kit (Ready to Use) (orb554510)

Description	Rabbit Matrix Metalloproteinase 3 (MMP3) Sandwich ELISA KitReady to Use: Pre-diluted detection reage...
Reactivity	Rabbit
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P28863
Sensitivity	0.061ng/mL
Entrez	100009111
Expiration Date	Please enquire.