

Product Datasheet

Rabbit MMP9 ELISA Kit (Ready to Use) (orb554511)

Description	Rabbit Matrix Metalloproteinase 9 (MMP9) Sandwich ELISA KitReady to Use: Pre-diluted detection reage...
Reactivity	Rabbit
Range	3.12-200ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P41246
Sensitivity	1.18ng/mL
Entrez	100008993
Expiration Date	Please enquire.