

Product Datasheet

Rabbit TIMP1 ELISA Kit (Ready to Use) (orb554533)

Description	Rabbit Tissue Inhibitors Of Metalloproteinase 1 (TIMP1) Sandwich ELISA KitReady to Use: Pre-diluted ...
Reactivity	Rabbit
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P20614
Sensitivity	0.119ng/mL
Entrez	100009047
Expiration Date	Please enquire.