

Product Datasheet

Rabbit TNNT2 ELISA Kit (Ready to Use) (orb554538)

Description	Rabbit Troponin T Type 2, Cardiac (TNNT2) Sandwich ELISA KitReady to Use: Pre-diluted detection reagent...
Reactivity	Rabbit
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P09741
Sensitivity	6.2pg/mL
Entrez	100009428
Expiration Date	Please enquire.