

Product Datasheet

Rat ADAM17 ELISA Kit (Ready to Use) (orb554551)

Description	Rat A Disintegrin And Metalloprotease 17 (ADAM17) Sandwich ELISA KitReady to Use: Pre-diluted detect...
Reactivity	Rat
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q9Z1K9
Sensitivity	0.112ng/mL
Entrez	57027
Expiration Date	Please enquire.