

Product Datasheet

Rat ADAMTS5 ELISA Kit (Ready to Use) (orb554557)

Description	Rat A Disintegrin And Metalloproteinase With Thrombospondin 5 (ADAMTS5) Sandwich ELISA KitReady to U...
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q6TY19
Sensitivity	0.064ng/mL
Entrez	304135
Expiration Date	Please enquire.