

Product Datasheet

Rat ACAT1 ELISA Kit (orb554559)

Description	Rat Acetyl Coenzyme A Acetyltransferase 1 (ACAT1) Sandwich ELISA Kit
Reactivity	Rat
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Sensitivity	12.5pg/mL
Expiration Date	Please enquire.