

Product Datasheet

Rat ACTb ELISA Kit (Ready to Use) (orb554567)

Description	Rat Actin Beta (ACTb) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for convenient ...
Reactivity	Rat
Range	1.56-100ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P60711
Sensitivity	0.64ng/mL
Entrez	81822
Expiration Date	Please enquire.