

Product Datasheet

Rat ADIPOR2 ELISA Kit (Ready to Use) (orb554575)

Description	Rat Adiponectin Receptor 2 (ADIPOR2) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents ...
Reactivity	Rat
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q3YAF7
Sensitivity	0.125ng/mL
Entrez	312670
Expiration Date	Please enquire.