

Product Datasheet

Rat ADH1 ELISA Kit (Ready to Use) (orb554588)

Description	Rat Alcohol Dehydrogenase 1 (ADH1) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents fo...
Reactivity	Rat
Range	1.56-100ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P06757
Sensitivity	0.59ng/mL
Entrez	24172
Expiration Date	Please enquire.